

ANALYTICAL SERVICES

A P C C A C O M P A N Y

LABORATORY REPORT

4/27/2015

Client #:

Sample: Pentobarbital Na5
50mg/ml

Lot #:

Sample ID #:

Date Rec'd: 4/14/2015

LABORATORY TEST RESULTS**Microbiological Tests:**

Bacterial Endotoxin USP <85>

<u>Date</u>	<u>Measured</u>	<u>Result</u>
4/17/2015	<1 EU/mL	Pass

Eagle Sterility Test

Day 7 (Neg.)

Rapid ScanRDI Microbial Detection

Chemical Tests:

Pentobarbital Sodium

<u>Date</u>	<u>Reported</u>	<u>Measured</u>	<u>Potency</u>
4/17/2015	50 mg/mL	49.9 mg/mL	99.8 %

Notes:


Endotoxins are measured using USP<85> Turbidimetric procedure, with an inhibition / enhancement test performed on each sample.

Eagle Sterility Membrane Filtration is a 14 day test using the membrane filtration procedure. Test for BACTERIA, MOLDS, YEAST AND FUNGI with two media at separate temperatures. Test read on days 3, and 7 are not final until the full 14 day test is complete.

Potency is determined via USP <621> HPLC, USP<851> Spectrophotometry, and specific monograph testing procedures.

Respectfully submitted,

EAGLE ANALYTICAL SERVICES LTD.



William J. Zolner, Ph.D., Chief Scientific Officer

9940 W. Sam Houston Pkwy S., Suite 310, Houston, TX 77099

832-295-1276

281-754-4700 (fax)

5000

EXHIBIT 3